

Product Datasheet

Recombinant Human Fc gamma RI / CD64 Protein, Avi His Tag (orb1550527)

Description	Recombinant Human Fc gamma RI / CD64 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	RI
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP1962, Recombinant Human Fc gamma RI / CD64 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gln16-Pro288) of Human Fc gamma RI / CD64 fused with a His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.