

Product Datasheet

Recombinant Human FGF10 Protein (orb1550528)

Description	Recombinant Human FGF10 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	FGF10
Note	For research use only
Purity	> 90% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1961, Recombinant Human FGF10 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Gln38-Ser208) of Human FGF10 fused with no tag .
Expiration Date	6 months from date of receipt.