

Product Datasheet

Recombinant Human EGFR VIII Protein, Avi His Tag (orb1550531)

Description	Recombinant Human EGFR VIII Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	EGFR VIII
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1958, Recombinant Human EGFR VIII Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Leu25-Ser645) of Human EGFR VIII fused with His Tag and Avi tag at the N-terminal.
Expiration Date	6 months from date of receipt.