

Product Datasheet

Recombinant Human FGF21 Protein, mFAvi (orb1550533)

Description	Recombinant Human FGF21 Protein, mFAvi
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	FGF21
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1956, Recombinant Human FGF21 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (His Tag29-Ser210) of Human FGF21 fused with a mFc and Avi tag at the N-terminal.
Expiration Date	6 months from date of receipt.