

Product Datasheet

Recombinant Human EGFR/HER1 Protein, hFc His Tag (orb1550536)

Description	Recombinant Human EGFR/HER1 Protein, hFc His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	EGFR
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1953, Recombinant Human EGFR/HER1 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Leu25-Ser645) of Human EGFR/HER1 fused with hFc His Tag tag at the C-terminal.
Expiration Date	6 months from date of receipt.