

Product Datasheet

Biotinylated Recombinant Human E-selectin Protein, Avi His Tag (orb1550537)

Description	Biotinylated Recombinant Human E-selectin Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Eselectin
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1952, Biotinylated Recombinant Human E-selectin Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Trp22-Pro556) of Human E-selectin fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.