

Product Datasheet

Biotinylated Recombinant Human EGFR VIII Protein, Avi His Tag (orb1550540)

Description	Biotinylated Recombinant Human EGFR VIII Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	EGFR VIII
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP1949, Biotinylated Recombinant Human EGFR VIII Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Leu25-Ser645) of Human EGFR VIII fused with His Tag and Avi tag at the N-terminal.
Expiration Date	6 months from date of receipt.