

Product Datasheet

Recombinant Human EGF Protein (orb1550544)

Description	Recombinant Human EGF Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	EGF
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1945, Recombinant Human EGF Protein is produced by E.coli expression system. The target protein is expressed with sequence (Asn971-Arg1023) of Human EGF fused with no tag .
Expiration Date	6 months from date of receipt.