

## Product Datasheet

# Recombinant Human CD24 Protein, Fc His Tag (orb1550549)

<b>Description</b>	Recombinant Human CD24 Protein, Fc His Tag
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	CD24
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
<b>Source</b>	KMP1940, Recombinant Human CD24 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ser27-Gly59) of Human CD24 fused with Llama IgG2b Fc tag at the C- terminal.
<b>Expiration Date</b>	6 months from date of receipt.