

Product Datasheet

Recombinant Cynomolgus CD3 Epsilon Protein, hFc His Tag (orb1550560)

Description	Recombinant Cynomolgus CD3 Epsilon Protein, hFc His Tag
Reactivity	Monkey
Endotoxins	< 1 EU/ug by LAL method.
Target	CD3
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP1929, Recombinant Cynomolgus CD3 Epsilon Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gln22-Asp117) of Cynomolgus CD3 Epsilon fused with a hFc His Tag tag at the C-terminal.
Expiration Date	6 months from date of receipt.