

Product Datasheet

Biotinylated Recombinant Human CD5 Protein, hFAvi His Tag (orb1550562)

Description	Biotinylated Recombinant Human CD5 Protein, hFAvi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CD5
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1927, Biotinylated Recombinant Human CD5 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Met1-Asn371) of Human CD5 fused with hFc tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.