

## Product Datasheet

# Recombinant Human CD52 Protein, mFc His Tag (orb1550565)

<b>Description</b>	Recombinant Human CD52 Protein, mFc His Tag
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	CD52
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
<b>Source</b>	KMP1924, Recombinant Human CD52 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gly25-Ser36 ) of Human CD52 fused with hFc His Tag tag at the C-terminal.
<b>Expiration Date</b>	6 months from date of receipt.