

Product Datasheet

Recombinant Human BCMA/TNFRSF17 Protein, hFc His Tag (orb1550567)

Description	Recombinant Human BCMA/TNFRSF17 Protein, hFc His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	BCMA
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1922, Recombinant Human BCMA/TNFRSF17 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Met1-Ala54) of Human BCMA/TNFRSF17 fused with a hFc His Tag tag at the C-terminal.
Expiration Date	6 months from date of receipt.