

## Product Datasheet

# Recombinant Human CD23 / Fc epsilon RII Protein, His Tag (orb1550572)

<b>Description</b>	Recombinant Human CD23 / Fc epsilon RII Protein, His Tag
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	CD23
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
<b>Source</b>	KMP1917, Recombinant Human CD23 / Fc epsilon RII Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Asp48-Ser321) of Human CD23 / Fc epsilon RII fused with a His Tag at the N-terminal.
<b>Expiration Date</b>	6 months from date of receipt.