

Product Datasheet

Recombinant Human CD122/IL2RB Protein, Avi His Tag (orb1550580)

Description	Recombinant Human CD122/IL2RB Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CD122
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1909, Recombinant Human CD122/IL2RB Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ala 27-Asp 239) of Human CD122/IL2RB fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.