

Product Datasheet

Recombinant Human Biotinyalted BCMA/TNFRSF17 Protein, Avi His Tag (orb1550582)

Description	Recombinant Human Biotinyalted BCMA/TNFRSF17 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	BCMA
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1907, Recombinant Human Biotinyalted BCMA/TNFRSF17 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Met1-Ala54) of Human Biotinyalted BCMA/TNFRSF17 fused with a His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.