

Product Datasheet

Recombinant Cynomolgus BAFFR Protein, His Tag (orb1550586)

Description	Recombinant Cynomolgus BAFFR Protein, His Tag
Reactivity	Monkey
Endotoxins	< 1 EU/ug by LAL method.
Target	BAFFR
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP1903, Recombinant Cynomolgus BAFFR Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ser7-Ala71) of Cynomolgus BAFFR fused with a His Tag at the C-terminal.
Expiration Date	6 months from date of receipt.