

Product Datasheet

Recombinant Human CA9 Protein, Avi His Tag (orb1550587)

Description	Recombinant Human CA9 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CA9
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP1902, Recombinant Human CA9 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gln38-Asp414) of Human CA9 fused with His Tag and Avi tag at the C- terminal.
Expiration Date	6 months from date of receipt.