

Product Datasheet

Recombinant Human PD-L1 Protein, mFc His Tag (orb1550594)

Description	Recombinant Human PD-L1 Protein, mFc His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	PDL1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1895, Recombinant Human CD30 Ligand/TNFSF8 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Phe19-Arg238) of human PD-L1 (Accession #Q9NZQ7) fused with an Fc tag at the C-terminus.
Expiration Date	6 months from date of receipt.