

## Product Datasheet

# Recombinant Mouse IL-1alpha/IL1A/IL1F1 Protein, His Tag (orb1550605)

<b>Description</b>	Recombinant Mouse IL-1alpha/IL1A/IL1F1 Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	IL1alpha
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1884, Recombinant Mouse IL-1alpha/IL1A/IL1F1 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ser115-Ser270) of Mouse IL1a (Accession #NP_034684.2) fused with an 6xHis Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.