

Product Datasheet

Recombinant Rat VEGFA/VEGF164 Protein, His Tag (orb1550607)

Description	Recombinant Rat VEGFA/VEGF164 Protein, His Tag
Species/Host	Human
Reactivity	Rat
Endotoxins	< 0.1 EU/ug by LAL method.
Target	VEGFA
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1882, Recombinant Rat VEGFA/VEGF164 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ala27-Arg190) of Rat VEGFA (Accession #NP_001274037.1) fused with an 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.