

Product Datasheet

Recombinant Human GDNF Protein with His tag, His Tag (orb1550617)

Description	Recombinant Human GDNF Protein with His tag, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	GDNF
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1872, Recombinant Human GDNF Protein with His Tag is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Arg109-Ile211) of human GDNF (Accession #NP_000505.1) fused with an 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.