

## Product Datasheet

# Recombinant MERS-CoV Spike S1 Protein with His tag, His Tag (orb1550618)

<b>Description</b>	Recombinant MERS-CoV Spike S1 Protein with His tag, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Virus
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	S1
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1871, Recombinant MERS-CoV Spike S1 Protein with His Tag is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Tyr18-Glu725) of mers-cov Spike S1 (Accession #AFS88936.1) fused with an 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.