

## Product Datasheet

# Recombinant Cynomolgus Cynomolgus CD3 Epsilon Protein, His Tag (orb1550631)

<b>Description</b>	Recombinant Cynomolgus Cynomolgus CD3 Epsilon Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Monkey
<b>Endotoxins</b>	< 1 EU per ug by the LAL method.
<b>Target</b>	CD3
<b>Note</b>	For research use only
<b>Purity</b>	> 95%by Bis-Tris PAGE, > 95%by HPLC
<b>Source</b>	KMP1858, Recombinant Cynomolgus Cynomolgus CD3 Epsilon Protein is produced by mammalian expression system. The target protein is expressed with sequence (Gln22-Asp117) of cynomolgus Cynomolgus CD3 Epsilon (Accession #) fused with a His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.