

Product Datasheet

Recombinant Human TNFSF13B Protein, Fc His Tag (orb1550633)

Description	Recombinant Human TNFSF13B Protein, Fc His Tag
Species/Host	Human
Reactivity	Human
Target	TNFSF13B
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1856, Recombinant Human TNFSF13B Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (hFc(Glu99-Ala330)&BAFF(Ala134-Leu285)) of human TNFSF13B (Accession #) fused with a Fc Tag at the N-terminal.
Expiration Date	6 months from date of receipt.