

## Product Datasheet

# Recombinant Human SDC1 Protein, Fc & His Tag (orb1550634)

<b>Description</b>	Recombinant Human SDC1 Protein, Fc & His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Target</b>	SDC1
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1855, Recombinant Human SDC1 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (CD138(Gln23-Gly254)&hFc(Glu99-Ala330)) of human SDC1 (Accession #) fused with a Fc and His Tag at the C-terminal.
<b>Expiration Date</b>	6 months from date of receipt.