

Product Datasheet

Recombinant Mouse CD3E&CD3D Protein, hFc His Tag (orb1550635)

Description	Recombinant Mouse CD3E&CD3D Protein, hFc His Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 0.1 EU per ug by the LAL method.
Target	CD3E
Note	For research use only
Purity	> 95%by Bis-Tris PAGE, > 95%by HPLC
Source	KMP1854, Recombinant Mouse CD3E&CD3D Protein is produced by mammalian expression system. The target protein is expressed with sequence (Asp23-Asp108(P22646)&Phe22-Ala105(P04235)) of mouse CD3E&CD3D (Accession #) fused with a hFc tag at the C-terminus.
Expiration Date	6 months from date of receipt.