

Product Datasheet

Recombinant Human TNFSF9 Protein, Fc & His Tag (orb1550636)

Description	Recombinant Human TNFSF9 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Target	TNFSF9
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1853, Recombinant Human TNFSF9 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (mFc(Pro99-Lys330)&4-1BB Ligand(Pro52-Glu254)) of human TNFSF9 (Accession #) fused with a Fc and His Tag at the N-terminal.
Expiration Date	6 months from date of receipt.