

Product Datasheet

Recombinant Human CD3 epsilon & CD3 delta Heterodimer Protein, hFAvi His Tag (orb1550643)

Description	Recombinant Human CD3 epsilon & CD3 delta Heterodimer Protein, hFAvi His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.05 EU per ug by the LAL method.
Target	CD3
Note	For research use only
Purity	> 95%by Bis-Tris PAGE, > 95%by HPLC
Source	KMP1846, Recombinant Human CD3 epsilon & CD3 delta Heterodimer Protein is produced by mammalian expression system. The target protein is expressed with sequence (Asp23-Asp126(CD3E) & Phe22-Ala105(CD3D)) of human CD3 epsilon & CD3 delta Heterodimer (Accession #) fused with a Fc and Avi Tag at the C-terminal.
Expiration Date	6 months from date of receipt.