

Product Datasheet

Recombinant Human GFRAL Protein, Avi His Tag (orb1550653)

Description	Recombinant Human GFRAL Protein, Avi His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1 EU per ug by the LAL method.
Target	GFRAL
Note	For research use only
Purity	> 95%by Tris-Bis PAGE, > 95%by HPLC
Source	KMP1836, Recombinant Human GFRAL Protein is produced by mammalian expression system. The target protein is expressed with sequence (Ser19-Glu351) of human GFRAL (Accession #) fused with a His Tag and Avi Tag at the C-terminal.
Expiration Date	6 months from date of receipt.