

Product Datasheet

Recombinant Human HAVCR2 Protein, Fc & His Tag (orb1550655)

Description	Recombinant Human HAVCR2 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Target	HAVCR2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1834, Recombinant Human HAVCR2 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (TIM3(Ser22-Gly202)&mFc(Pro99-Lys330)) of human HAVCR2 (Accession #) fused with a Fc and His Tag at the C-terminal.
Expiration Date	6 months from date of receipt.