

Product Datasheet

Recombinant Human TNFSF8 Protein, Fc & His Tag (orb1550658)

Description	Recombinant Human TNFSF8 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Target	TNFSF8
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1831, Recombinant Human TNFSF8 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (mFc(Pro99-Lys330)&CD30 Ligand(Gln63-Asp234)) of human TNFSF8 (Accession #) fused with a Fc and His Tag at the N-terminal.
Expiration Date	6 months from date of receipt.