

## Product Datasheet

# Recombinant Human Siglec10 Protein, His Tag (orb1550664)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Recombinant Human Siglec10 Protein, His Tag   |
| <b>Species/Host</b>    | Human   |
| <b>Reactivity</b>      | Human   |
| <b>Target</b>          | Siglec10  |
| <b>Note</b>            | For research use only   |
| <b>Purity</b>          | > 95% by SDS-PAGE.  |
| <b>Source</b>          | KMP1825, Recombinant Human Siglec10 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (SIGLEC10(Met17-Thr546)) of human Siglec10 (Accession #) fused with a His Tag at the C-terminus. |
| <b>Expiration Date</b> | 6 months from date of receipt.  |