

Product Datasheet

Recombinant Human SIGLEC15 Protein, Fc & His Tag (orb1550669)

Description	Recombinant Human SIGLEC15 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Target	SIGLEC15
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1820, Recombinant Human SIGLEC15 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (SIGLEC15(Phe20-Thr263)&mFc(Pro99-Lys330)) of human SIGLEC15 (Accession #) fused with a Fc and His Tag at the C-terminal.
Expiration Date	6 months from date of receipt.