

Product Datasheet

Recombinant Human IL3RA Protein, Fc & His Tag (orb1550670)

Description	Recombinant Human IL3RA Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Target	IL3RA
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1819, Recombinant Human IL3RA Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (CD123(Thr19-Arg305)&hFc(Glu99-Ala330)) of human IL3RA (Accession #) fused with a Fc and His Tag at the C-terminal.
Expiration Date	6 months from date of receipt.