

## Product Datasheet

# Recombinant Human IL2RB Protein, His Tag (orb1550671)

<b>Description</b>	Recombinant Human IL2RB Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Target</b>	IL2RB
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1818, Recombinant Human IL2RB Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (IL2RB(Ala27-Asp239)) of human IL2RB (Accession #) fused with a His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.