

Product Datasheet

Recombinant Human CD37 Protein, Fc His Tag (orb1550675)

Description	Recombinant Human CD37 Protein, Fc His Tag
Species/Host	Human
Reactivity	Human
Target	CD37
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1814, Recombinant Human CD37 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (CD37(Arg112-Asn241)&hFc(Glu99-Ala330)) of human CD37 (Accession #) fused with a Fc Tag at the C-terminal.
Expiration Date	6 months from date of receipt.