

Product Datasheet

Recombinant Human SELPLG Protein, His Tag (orb1550678)

Description	Recombinant Human SELPLG Protein, His Tag
Species/Host	Human
Reactivity	Human
Target	SELPLG
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1811, Recombinant Human SELPLG Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (CD162(Leu18-Cys320)) of human SELPLG (Accession #) fused with a His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.