

Product Datasheet

Recombinant Human IL5RA Protein, His Tag (orb1550680)

Description	Recombinant Human IL5RA Protein, His Tag
Species/Host	Human
Reactivity	Human
Target	IL5RA
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1809, Recombinant Human IL5RA Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (IL5RA(Asp21-Pro338)) of human IL5RA (Accession #) fused with a His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.