

Product Datasheet

Recombinant Human CD34 Protein, His Tag (orb1550681)

Description	Recombinant Human CD34 Protein, His Tag
Species/Host	Human
Reactivity	Human
Target	CD34
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1808, Recombinant Human CD34 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (CD34(Ser32-Thr290)) of human CD34 (Accession #) fused with a His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.