

## Product Datasheet

# Recombinant MERS-CoV Nucleocapsid Protein with His tag, His Tag (orb1550685)

<b>Description</b>	Recombinant MERS-CoV Nucleocapsid Protein with His tag, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Virus
<b>Endotoxins</b>	Please contact us for more information.
<b>Target</b>	MERSCoV
<b>Note</b>	For research use only
<b>Purity</b>	> 90% by SDS-PAGE.
<b>Source</b>	KMP1804, Recombinant MERS-CoV Nucleocapsid Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Asp413) of mers-cov Nucleocapsid (Accession #YP_009047211.1) fused with a 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.