

## Product Datasheet

# Recombinant Mouse Siglec-15/CD33L3 Protein with Fc tag, Fc His Tag (orb1550692)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Recombinant Mouse Siglec-15/CD33L3 Protein with Fc tag, Fc His Tag   |
| <b>Species/Host</b>    | Human  |
| <b>Reactivity</b>      | Mouse  |
| <b>Endotoxins</b>      | < 0.1 EU/ug by LAL method.   |
| <b>Target</b>          | Siglec15   |
| <b>Note</b>            | For research use only  |
| <b>Purity</b>          | > 90% by SDS-PAGE.   |
| <b>Source</b>          | KMP1797, Recombinant Mouse Siglec-15/CD33L3 Protein with Fc tag is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Arg24-Thr262) of mouse Siglec-15/CD33L3 (Accession #NP_001094508.1) fused with a Fc tag at the C-terminus. |
| <b>Expiration Date</b> | 6 months from date of receipt.   |