

Product Datasheet

Recombinant Mouse Siglec-15/CD33L3 Protein with His tag, His Tag (orb1550694)

| | |
|------------------------|--|
| Description | Recombinant Mouse Siglec-15/CD33L3 Protein with His tag, His Tag |
| Species/Host | Human |
| Reactivity | Mouse |
| Endotoxins | < 1.0 EU/ug by LAL method. |
| Target | Siglec15 |
| Note | For research use only |
| Purity | > 90% by SDS-PAGE. |
| Source | KMP1795, Recombinant Mouse Siglec-15/CD33L3 Protein with His Tag is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Arg24-Thr262) of mouse Siglec-15/CD33L3 (Accession #NP_001094508.1) fused with a 6×His Tag at the C-terminus. |
| Expiration Date | 6 months from date of receipt. |