

## Product Datasheet

# Recombinant Human ACE2 Protein with His tag, His Tag (orb1550695)

<b>Description</b>	Recombinant Human ACE2 Protein with His tag, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	ACE2
<b>Note</b>	For research use only
<b>Purity</b>	> 90% by SDS-PAGE.
<b>Source</b>	KMP1794, Recombinant Human ACE2 Protein with His Tag is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gln18-Ser740) of human ACE2 (Accession #Q9BYF1) fused with a 6×His Tag at the N-terminus.
<b>Expiration Date</b>	6 months from date of receipt.