

Product Datasheet

Recombinant Human ACE2 Protein with hFc tag, Fc His Tag (orb1550698)

Description	Recombinant Human ACE2 Protein with hFc tag, Fc His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	ACE2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1791, Recombinant Human ACE2 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gln18-Ser740) of human ACE2 (Accession #Q9BYF1) fused with a Fc tag at the C-terminus.
Expiration Date	6 months from date of receipt.