

Product Datasheet

Recombinant human S100A9 Protein, His Tag (orb1550700)

Description	Recombinant human S100A9 Protein, His Tag
Species/Host	Virus
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	S100A9
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1789, Recombinant human S100A9 Protein is produced by Baculovirus-Infected Sf9 Cells expression system. The target protein is expressed with sequence (Met1-Pro114) of human S100A9 (Accession #NP_002956.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.