

Product Datasheet

Recombinant human IL-15 Protein, Fc & His Tag (orb1550702)

Description	Recombinant human IL-15 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	IL15
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1787, Recombinant human IL-15 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Asn49-Ser162) of human IL-15 (Accession #NP_000576.1) fused with a Fc, 6xHis Tag at the C- terminus.
Expiration Date	6 months from date of receipt.