

Product Datasheet

Recombinant human S100A8/S100A9 Heterodimer Protein, His Tag (orb1550703)

Description	Recombinant human S100A8/S100A9 Heterodimer Protein, His Tag
Species/Host	Virus
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	S100A8
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1786, Recombinant human S100A8/S100A9 Heterodimer Protein is produced by Baculovirus-Infected Sf9 Cells expression system. The target protein is expressed with sequence (Met1-Glu93&Met1-Pro114) of human S100A8&S100A9 Heterodimer (Accession #NP_002955.2/NP_002956.1) fused with a 6xHis Tag at the C-terminus(S100A8).
Expiration Date	6 months from date of receipt.