

Product Datasheet

Recombinant Human IL-7 Protein, His Tag (orb1550704)

Description	Recombinant Human IL-7 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	IL7
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1785, Recombinant Human IL-7 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Asp26-His Tag177) of human IL-7 (Accession #NP_000871.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.