

Product Datasheet

Recombinant Human S100 Calcium Binding Protein B/S100B Protein, His Tag (orb1550705)

Description	Recombinant Human S100 Calcium Binding Protein B/S100B Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	S100 Calcium Binding B
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1784, Recombinant Human S100 Calcium Binding Protein B/S100B Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Glu92) of human S100 Calcium Binding Protein B/S100B (Accession #P04271) fused with an initial Met at the N-terminus and a 6×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.